

Bauer Core Standard Protocol		
Title: Using the Covaris S220 to shear DNA in a micro tube		
Pages: 2	Revision: 2.0	Date: 3/9/16
Author(s): Claire Reardon	Reviewers:	
Contact: claire@cgr.harvard.edu	Comment:	

1. Purpose

This protocol provides instructions for using the Covaris S220 shearing instrument to shear nucleic acid samples. The protocol is designed as a reference and is not a substitute for training. Users must complete a training session before using any of the Bauer Core's instrumentation.

2. Materials

- 2.1. Sample(s) to be sheared
 - 2.1.1. Sample quantity: up to 10 µg
 - 2.1.2. Sample volume: 130µl or 50 µl
- 2.2. Microtube (6 x 16 mm) part # 520045

3. Instrumentation Covaris S220

4. Reagent preparation None.

5. Procedure

- 5.1. Start Up: Allow 30 minutes for start up procedure.
 - 5.1.1. Turn on the chiller beneath the bench (power switch on back).
 - 5.1.1.1. Set to 2.5 °C to give water at ~6°C in the Covaris water bath.
 - 5.1.1.2. To adjust, push the knob, turn to desired temp, and push again.
 - 5.1.2. Fill the Covaris water bath with DI H₂O.
 - 5.1.3. Run the SonoLab software.
 - 5.1.4. Allow the water bath to degas and come to temperature before running.
- 5.2. Creating a New Method
 - 5.2.1. Refer to library prep protocol for settings, or select according to the tables below:

For 130ul

Target BP	150	200	300	400	500	800	1000	1500
Duty Factor	10%	10%	10%	10%	5%	5%	5%	2%
Peak Incident Power	175	175	140	140	105	105	105	140
Cycles Per Burst	200	200	200	200	200	200	200	200
Duration	430	180	80	55	80	50	40	15

For 50ul

Target BP	150	200	300	400	500	1000	1500
Duty Factor	10%	10%	10%	5%	5%	2%	1%
Peak Incident Power	175	175	175	175	175	175	175
Cycles Per Burst	200	200	200	200	200	200	200
Duration	280	120	50	55	35	45	20

5.3. Preparing the Sample Tube

5.3.1. Load the sample through the septum into the tube with a pipette.

5.3.2. Tubes are designed for 130µl but can also accommodate 50 µl.

5.3.3. Check for air pockets and spin down if present.

5.3.4. Place the tube into the plunger holder so that it sits straight.

5.3.5. Load the plunger holder into the instrument (any direction).

5.4. Running a Saved Method.

5.4.1. Wait for the instrument status to be ready and then hit Run.

5.5. Cleanup

5.5.1. Remove the sample tube.

5.5.2. Turn off the degas pump

5.5.3. Empty the water tank.

5.5.4. Run the degas pump to purge residual water from the lines.

The pump will automatically stop after 10 seconds.

5.5.5. Empty residual water from the tank and allow to air dry.