

Bauer Core Standard Protocol		
Title: Using the Covaris S220		
Pages: 2	Revision: 1.1	Date: 1/26/2017
Author(s): Claire Reardon	Reviewers: Nicole El-Ali	
Contact: claire@cgr.harvard.edu	Comment: General protocol. See specific protocol for DNA shearing instructions.	

1. Purpose

To use the Covaris S220 shearing instrument. Users must complete a training session before using any Bauer Core instrumentation.

2. Materials

2.1. Sample(s) to be sheared

2.2. Tube for desired protocol

2.2.1. MicroTube: 150 - 1500 bp, 130µl or 50µl - part #520045 (provided by core)

2.2.2. MiniTube Clear: 2000 bp, 200ul - part # 520064 (provided by core)

2.2.3. MiniTube Blue: 3000 bp, 200ul - part # 520065 (provided by core)

2.2.4. MiniTube Red: 5000 bp, 200ul - part # 520066 (provided by core)

2.2.5. MilliTube: 1 ml - part # 520135 (NOT provided by core)

2.2.6. MicroTube-15: 150 - 1500bp, 15ul - part # 520145 (NOT provided by core)

3. Instrumentation

3.1. Covaris S220

3.2. Tube holder (tube holders available in core are listed below)

3.2.1. Microtube holder part # 500114

3.2.2. MiniTube holder part # 500206

3.2.3. MilliTube holder part # 500371

3.2.4. MicroTube-15 holder part # 500427

3.3. Microcentrifuge (minifuge may be used for MicroTubes, but NOT MicroTube-15)

3.3.1. Microcentrifuge adapter #500406

4. Reagent preparation None.

5. Procedure

5.1. Start Up: Allow 30 minutes for start up procedure.

5.1.1. Turn on the chiller beneath the bench (power switch on back).

5.1.1.1. Set chiller to desired temp. (2.5 °C gives ~7°C in Covaris tank).

5.1.1.2. To adjust, push the knob, turn to desired temp, and push again.

5.1.2. Fill the Covaris tank to the level desired with DI H₂O.

5.1.3. Turn on the Covaris instrument (power switch on front)..

- 5.1.4. Open the SonoLab software.
- 5.1.4. Allow the water bath to degas and come to temperature.
- 5.2. Create a method with desired parameters
 - 5.2.1. Standard DNA shearing methods are in the “Bauer Core” folder.
- 5.3. Preparing the Sample Tube
 - 5.3.1. Pipette the sample into the tube (many tubes have septa for easy loading).
 - 5.3.2. Check for air pockets and spin down if present.
 - 5.3.3. Place the tube into the holder so that it sits straight.
 - 5.3.4. Load the holder into the instrument (any direction).
- 5.4. Running a Saved Method
 - 5.4.1. Select a method (standard DNA shearing methods are in the “Bauer Core” folder).
 - 5.4.2. Wait for the instrument status to be ready and then hit Run. Wait for the instrument status to be ready and then hit Run.
- 5.5. Cleanup
 - 5.5.1. Remove the sample tube.
 - 5.5.1.1. Spin tube briefly to collect sample.
 - 5.5.1.2. Remove cap to pipette out sample.
 - 5.5.2. Turn off the degas pump
 - 5.5.3. Empty the water tank.
 - 5.5.4. Run the degas pump to purge residual water from the lines.
 - The pump will automatically stop after 10 seconds.
 - 5.5.5. Empty residual water from the tank and allow to air dry.