

**Before your run:**

- Purchase tape(s) [through minilims at least 24-48 hours before](#) your reservation time
  - Reference tape specs on the front of the instruments or [see here for an online version.](#)
- Reserve instrument in [SPINAL](#)
- Warm tape(s) and reagents to room temperature before starting instrument set-up (~30 mins)

**Instrument Set-Up:**

- Instrument, Reagents, Consumables:
  - Available on bench for users to use: pipettes, cyclers, microfuges, other items in the surrounding area
  - Plates & foil seals, strip tubes & caps, tips in drawer
    - **Always use** Agilent designated consumables, especially foil seal.
    - **Always remove caps from strip tubes prior to starting instrument**
    - **Do Not cut strip tubes**
  - Quick Guides are in the drawer -always use for reference
    - Volume of buffer & ladder/sample differs depending on tape used
    - See last pages for images of each set-up
  - Using the instrument
    - Power button & power cycling
      - Turn instrument on & then software, do reverse to turn off
    - Instrument deck set-up
      - Important to make sure tape waste & tip waste are always cleaned prior to starting run
  - Reagents located in mini-fridge below the middle bench
    - DNA has a ladder provided by Core
    - RNA does not
      - **Ladder will always be in position A1 of strip tube**
      - Ladder is run on every new tape for DNA - 15 samples per tape
      - RNA will need to use electronic ladder - 16 samples per tape
    - Order tapes through minilims & pick up in drop off fridge
    - Orient tape in instrument with the barcode in the bottom right corner, facing the back so that you see the blue or green label
      - Ensure no bubbles in lanes before use
        - Lightly tap to dislodge any bubbles in the lane
    - Partial tapes can't be used on a different instrument
      - Use color-coding stickers and labeling to ensure you do not mix your tapes
        - **Orange #1**
        - **Red #2**
    - When samples are set up, vortex **using designated vortexer set at 2000 rpm** & spin down

### Laptop log-in:

- Log in with RC account
- Open controller software
  - First time use will need to change settings to **allow expired tapes**
  - Can select wells for both strip tubes or plate format
    - Recall: DNA ladder is always in A1 location of the strip tube
    - For RNA, right click and select “Electronic Ladder”
  - Import .txt file for descriptions
    - Create an excel .txt file with the “sample description” in second column
      - No need to add in any well locations, it will automatically be set by file
    - **Save as .txt** to RC account to import at instrument – see details below for how to access the location
  - Prefix is the file name – unique to your labeling and organizing process
    - Don’t need to include date if it is part of output under settings
  - Required box will show how many tips, tapes, & ladder/buffer need for run
    - Important to doublecheck this prior to starting run
  - When ready select start, then proceed, & **wait for run to start**
    - If you walk away too soon, it could error out, so wait for the status bar

### Analysis software

- First under File, Save As to RC folder
  - This will save the original data file
- De-select align to check alignment
  - If there are issues see troubleshooting guide
- For DNA, view the samples under the Region tab
  - Also select scale to sample
  - Add region(s)
    - Can either right click over electropherogram and select “add region”  
OR
    - Can go to region settings and set a standard region across the entire file
- For RNA, stay under Electropherogram tab
  - May need to right click and select the sub unit peaks if not appropriately called
- Scroll through each profile to see if generally everything looks good
  - No major tape or sample related issues
- Can re-assign, add, remove peaks while under Electropherogram tab
- Export file
  - RNA: File → Export → for RNA select sample table only
  - DNA: File → Export → for DNA select compact region table only
- Create report
  - Select/deselect any samples in the “samples” location
  - Deselect Peak table for DNA – region table will have the necessary information

**When the run is complete:**

- Save partial tape
  - Wrap up with original foil seal and tape
  - Write the latest date on the foil so the core does not toss it out during a tape clean out – note: anything 4 weeks or older will be discarded as well as anything not labeled
  - Unused, old tapes will be recovered if older than 4 weeks
- Clean up bench area
  - Put things away and wipe down with 70% EtOH
  - Log out of your windows account, DO NOT lock the profile

**Accessing and Saving to Your RC Account: Data Storage and Transfer Policy**

**Note:** Data storage on local hard drives is prohibited, and any files found on the hard drives may be deleted without warning.

- **Users with lab folders on the Research Computing Network Storage Space:**
  - Can connect directly to their folders on RC file systems (\\rcfs1, \\rcfs2, \\rcfs3, etc...) and should use these servers to store their data.
  - There are no time limits for these servers as they belong to the individual labs.
- **Users without Lab Folders:**
  - Can temporarily store data generated on the center's instruments in a folder on [\\rcstore02.rc.fas.harvard.edu\data](https://rcstore02.rc.fas.harvard.edu/data) you can access this through the file explorer
  - Within this area, users can create folders and store up to 2.5 GB of data.
  - If you don't have a lab folder, this will be mapped as your Z: drive in Windows Explorer.
    - These data will be accessible to download remotely for 90 days by either HTTPS or SFTP. See instructions below.
    - You will need to authenticate with your Active Directory (RC) username and password.
- **HTTPS:**
  - Point your web browser to <https://data.rc.fas.harvard.edu>
  - Log in using your Research Computing username and your password. Open your subfolder and click on files to open or save them. This will only allow you to open/save a single file at a time.
- **SFTP:**
  - Follow the instructions at <https://rc.fas.harvard.edu/resources/access-and-login/#Odyssey> access requires the OpenAuth tool for two factor authentication to establish your OpenAuth two-factor authentication tool.
  - Use an FTP client such as Filezilla (freeware). Under Site Manager, create a new site (Host: data.rc.fas.harvard.edu)
    - Protocol: SFTP Logon Type: Interactive
  - Under Advanced, Default remote directory: /data
  - Under Transfer Settings, check the box for "Limit number of simultaneous connections" and set the maximum number of connections to 1.
    - See the following page for more details: <https://rc.fas.harvard.edu/resources/documentation/transferring-data/sftp-file-transfer/>
  - Browse for your data in the folder in the lower right field.
  - Browse for the desired destination in the lower left field.
  - Drag files from the right hand field to the left hand field to copy them from the data folder to your computer.