

LSR Fortessa Special Order System Instrument Configuration				
Laser	Parameter Name	Long pass filter	bandpass filter	PMT
440nm	mCFP	n/a	483/32 488 Notch	1
440nm	Chromomycin A3	570	585/42	2

488nm	SSC	n/a	488/10	3
488nm	FITC	505	530/30	4
488nm	PerCp	670	685/35	5

561 nm	PE	n/a	582/15	6
561 nm	PE Texas Red	600	610/20	7
561 nm	Pe Cy 5	635	660/20	8
561 nm	Pe Cy 7	750	780/60	9

635 nm	APC	n/a	670/30	10
635 nm	APC Cy 7	750	780/60	11

*FSC detector is a Photo Diode

Note: All optical filter holders have a YELLOW colored label.